# Medtronic

DBS 3387  
Lead kit for deep brain  
stimulation 3389

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







! USA Rx only

**CE 0123**  
1998


## Explanation of symbols on product or package labeling

Refer to the appropriate product for symbols that apply.

**CE 0123** Conformité Européenne (European Conformity). This symbol means that the device fully complies with European Directive AIMD 90/385/EEC.

Medtronic and Medtronic logo are trademarks of Medtronic. Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.

# Table of contents

**Description 7**

**Package contents 7**

**Device specifications 7**

**Instructions for use 10**

Preparing for surgery 10

Placing a DBS lead 10

Stylet removal and lead stabilization 12

Capping the lead 16

Removing the lead cap 19

## **Information available for the Deep Brain Stimulation Therapies:**

The information for prescribers manual provides information about contraindications, warnings, precautions, adverse events, sterilization, and component disposal. The information for prescribers addendums contain indication-specific information. For customers in Japan, the appropriate package insert provides information about safety, contraindications, warnings, precautions, and adverse events.

The indications sheet provides information about indications and related information. For customers in Japan, the appropriate package insert provides information about indications.

The system eligibility and battery longevity manual describes programming considerations and provides battery longevity information to aid in the appropriate neurostimulator selection.

MRI guidelines provide information about any MRI conditions and MRI-specific contraindications, warnings, and precautions for MRI scans with the neurostimulation system.

Product manuals, such as programming guides, recharging guides, and implant manuals provide device descriptions, package contents, device specifications, product-specific warnings and precautions, and instructions for use.

**USA** The clinical summary provides information about the clinical study results for the neurostimulation system.

## Description

The Medtronic Models 3387 and 3389 DBS Leads are part of a neurostimulation system for deep brain stimulation.


## Package contents

- One Model 3387 or 3389 Lead
- Stylets:
  - straight (inserted in lead)
  - short
- Torque wrench
- Depth stop gauge (lead)
- Burr hole ring and cap
- Connector boot
- Tunneling tools: tunneling rod, tunneling tip, tube (passing straw)
- Lead cap
- Product literature
-  Warranty card
-  Registration form

## Device specifications

The lead has electrodes on the distal end; the proximal end fits into a four-conductor connector. A stylet has been inserted into the proximal end of the lead to aid in positioning the lead and performing intraoperative test stimulation.

- The Model 3387 DBS Lead features wide (1.5 mm) spacing between each of the four electrodes at the distal end.
- The Model 3389 DBS Lead features narrow (0.5 mm) spacing between each of the four electrodes at the distal end.
- The physician may choose the Model 3387 DBS Lead to provide electrodes spread over 10.5 mm or the Model 3389 DBS Lead to provide electrodes spread over 7.5 mm.

**Figure 1.** Electrode, lead contact, and stylet.

**Table 1. Device specifications<sup>a</sup> for lead Models 3387 and 3389**

<b>Description</b>	<b>Model 3387</b>	<b>Model 3389</b>
Connector	Quadripolar, in-line	Quadripolar, in-line
Shape	Straight	Straight
Conductor resistance <sup>b</sup>	<100 Ω	<100 Ω
Length	10 – 50 cm	10 – 50 cm
Diameter	1.27 mm	1.27 mm
Distal end		
Number of electrodes	4	4
Electrode shape	Cylindrical	Cylindrical
Electrode length	1.5 mm	1.5 mm
Electrode spacing	1.5 mm	0.5 mm
Electrode distance	10.5 mm	7.5 mm
Distal tip distance	1.5 mm	1.5 mm
Proximal end		
Lead contact length	2.3 mm	2.3 mm
Lead contact spacing	4.3 mm	4.3 mm
Lead contact distance	16.6 mm	16.6 mm
Stylet handle length	40.1 mm	40.1 mm

<sup>a</sup> All measurements are approximate.

<sup>b</sup> Electrical resistance is proportional to lead length.

**Table 2. Material of components in the Models 3387 or 3389 package, or in the Models 3550-03 or 3550-07 package**

<b>Component</b>	<b>Material</b>	<b>Material contacts human tissue</b>
Lead		
Conductor wire	Platinum-iridium	No
Conductor wire insulation	Fluoropolymer	Yes
Electrodes	Platinum-iridium	Yes
Stylet insulation	Parylene	Yes
Outer jacket tubing	Polyurethane	Yes
Proximal connector	MP35N	Yes
Lead cap		
Lead cap	Silicone rubber	Yes
Setscrew connector block	Stainless steel	Yes
Setscrew	Titanium	Yes
Connector boot	Silicone rubber	Yes
Burr hole ring and cap		
Ring	Polyamide	Yes
Cap	Silicone rubber	Yes
TwistLock screening cable (in the Model 3550-03 package)		

**Table 2. Material of components in the Models 3387 or 3389 package, or in the Models 3550-03 or 3550-07 package** (continued)

Component	Material	Material contacts human tissue
Connector	Polycarbonate	Yes
Cable	Rubber	Yes
Alligator clip screening cable (in the Model 3550-07 package)		
Connector	Polypropylene	Yes
Cable	Rubber	Yes
Clip	Stainless steel	Yes
Boot	Thermo plastic rubber (Santoprene™)	Yes
Stylets (straight, short)		
Handle	Acetal resin	Yes
Wire	Tungsten	Yes
Insulation	Fluoropolymer	Yes
Tunneling tools		
Tunneling rod	Stainless steel	Yes
Tunneling tip	Stainless steel	Yes
Tube (passing straw)	Polytetrafluoroethylene	Yes
Torque wrench		
Handle	Polypropylene	Yes
Shaft	Stainless steel	Yes

## Instructions for use

Implanting physicians should be experienced in stereotactic and functional neurosurgery and deep brain stimulation procedures, as well as thoroughly familiar with all product labeling.

The deep brain stimulation system delivers electrical stimulation to selected targets in the brain. The target site may be localized for stereotactic implantation of the DBS lead using CT scans, MRI, or ventriculography. A test stimulation or mapping electrode may then be utilized for further localization of the target. The DBS lead is not recommended for mapping.

### Preparing for surgery

Before opening the lead package, verify the model number, use-by date, lead length, and connector type.

### Placing a DBS lead

Medtronic recognizes that a variety of techniques may be used for the initial implant. The following procedure, which uses a stereotactic frame and an insertion cannula, is presented as one possible approach for the physician's consideration.

**Warning:** If damage occurs (eg, breaking, cutting, nicking, flattening, stretching, etc.) to the lead or any lead-accessory component during the procedure, do not implant the lead. Use a new lead.

### Cautions:

- The stylet inserted in the Model 3387 or 3389 Lead is specific to each lead; stylets are not interchangeable between lead models or between individual leads of the same model. Use of a nonmatching stylet may damage the lead, which would require surgical replacement.
- Do not bend, kink, stretch, or twist the lead or extension, which may damage the component. Component damage would require surgical replacement.
- Do not bend or kink the stylet, which may damage the component. Component damage would require surgical replacement.
- Ensure that the lead contacts are aligned with the stylet handle when connected to the test stimulator cable. If the lead contacts are not aligned properly, intermittent or loss of stimulation may occur.
- Use only rubber-tipped forceps on the lead. Do not use sharp-edged instruments (eg, hemostat), which may nick or cut the insulation, and require surgical replacement of the lead.
- Do not use saline or other ionic fluids at connections, which could result in a short circuit.

**Note:** During a lead implant procedure for a movement disorder, an insertion cannula should be placed to a point approximately 15 mm proximal to the target site for stimulation. The lead should be passed through the insertion cannula and advanced to the target site. A stereotactic frame kit may be used to stabilize the lead in the cannula.

1. After placement of the stereotactic frame, use standard imaging techniques to determine coordinates for the lead target site.
2. Prepare the patient per normal stereotactic neurosurgical techniques.
3. Make a skin incision, with consideration given to burr hole placement.
4. Prepare a subgaleal pocket by blunt dissection at the top of the skull at the edge of the burr hole incision for placement of the excess lead wire and connector.

**Note:** Placement of the pocket may be on either the left or right side of the skull.

5. Create a 14-mm diameter burr hole in the desired location.

### Notes:

- Medtronic recommends using a straight-edged 14-mm perforator to form the burr hole.
- Refer to the technical manual provided with the StimLoc™ burr hole cover for instructions for use on placing the burr hole cover.


**Warning:** Use only the burr hole cap and ring provided in the lead kit or Model 924256 StimLoc burr hole cover to anchor the lead. Do not use any other burr hole cover or alternative methods to anchor the lead. Use of any other anchoring method for example, a different model burr hole cover, tying ligatures, glues, cements, surgical plates, etc. may damage the lead, which would require surgical replacement.

6. Place the burr hole ring tightly against the bone in the burr hole, using your finger and a curved mosquito hemostat.
7. Attach the lead holder assembly from the frame kit to the stereotactic frame, positioning the guide tube or collimator so that its distal end is 1.25 to 2.5 cm from the skull (Figure 2).

**Note:** Refer to the stereotactic frame kit manual for instructions.

- a. Insert the stylet into the cannula.
- b. Advance the insertion cannula through the guide tube to a point approximately 15 mm proximal to the lead target site.
- c. Remove stylet from cannula.

**Figure 2.** Attach lead holder assembly to frame.

8. Determine target location with a test electrode.
  - a. Place the test electrode 10 mm above the target. Slowly advance it 1 mm at a time.
  - b. Once the target is determined, document the location and remove the test electrode.

**Note:** It may be necessary to make additional tracks to achieve optimum lead location. However, the risks of additional tracks should be weighed against potential benefits.

9. Determine the depth to place the lead, and set the delivery system to the proper target depth.
10. Measure the lead to the appropriate depth and attach the depth stop included in the lead kit.
11. Insert the lead through the insertion cannula to the measured target depth. Adjust the lead holder to less than 2.5 cm below the stylet handle, then attach the lead to the lead holder. The lead should be straight to prevent migration during lead placement.
12. Advance the lead slowly along the track made by the test electrode.

**Caution:** An increase in resistance or friction during lead insertion may indicate that the lead is deviating from the intended track. If this occurs, pull the lead back and re-advance until the target is reached.


13. Perform lead targeting per the appropriate indication for use.  
Maximum therapeutic effect with minimum side effects indicate the appropriate target.

## Stylet removal and lead stabilization

Complete the following steps to secure the positioned lead in the burr hole ring.


1. Ensure that the lead is securely fastened in the lead holder groove. The lead holder is provided with the lead frame kit.
  - a. Remove the adjustable depth stop gauge from the lead, taking care not to dislodge the lead.

- b. Carefully pull the insertion cannula up until the lead can be seen between the burr hole and the cannula.
- c. Gently hold the lead at the point it exits the skull. An assistant should hold the insertion cannula to ensure that it stays in position.
- d. Carefully loosen the stylet handle from the lead (Figure 3).
- e. While keeping the lead secure at the exit site, remove the lead from the lead holder.
- f. Remove the stylet from the lead.
- g. Remove the insertion cannula.
- h. Remove the guide tube assembly.

**Figure 3.** Loosen the stylet handle from the lead.


2. Gently press the lead into one of the precut grooves on the inner side of the burr hole ring, ensuring that the lead exits the burr hole ring on the side which the neurostimulator will be implanted (Figure 4).

**Figure 4.** Gently press the lead into the burr hole ring.


3. To minimize the potential for lead movement, place the burr hole cap into the ring as follows:
  - a. Align the tab in the burr hole cap with the slot in the burr hole ring.
  - b. Hold the cap tilted toward the slot in the ring where the lead is fixed, and gently press the cap against the edge of the burr hole ring and lead (Figure 5a).

- c. Gently press down on the lead exit side of the burr hole cap while rolling the cap downwards into the ring until secure (Figure 5b).

**Figure 5.** Gently roll burr hole cap into burr hole ring.


4. Recheck the stimulation effect after stabilizing the lead in the burr hole ring and cap. If using the alligator clip (packaged separately) do the following:
- Check that the test stimulator output (Amplitude) is adjusted to **off**.
  - Carefully attach the alligator clips to the desired connector ring contacts on the lead end.
5. If using the TwistLock screening cable (packaged separately) to recheck stimulation, complete the following steps:
- Attach the enclosed short stylet to the proximal end of the lead (Figure 6).

**Figure 6.** Attach short stylet to lead.


- Insert the short stylet completely into the lead (Figure 6a).
- Secure the lead in the stylet handle (Figure 6b).
- Insert the stylet handle of the lead into the TwistLock connector on the screening cable and lock it (refer to Figures 7 - 10).

**Note:** The handle fits into the TwistLock connector in only one way (Figure 7).


**Figure 7.** Insert the stylet handle into the TwistLock connector.


**Note:** Place the stylet handle into the groove at a slight angle to secure the handle end (Figure 8).

**Figure 8.** Place the stylet handle in the groove.

**Figure 9.** Press the stylet handle end in the groove.

**Figure 10.** Turn to lock the TwistLock connector.


6. Turn on the test stimulator and recheck stimulation effects.

**Note:** If lead movement has occurred, it may be necessary to remove the lead and repeat the implant procedure, using a new lead.


**Warning:** If lead repositioning is required, do not reinsert stylet into implanted lead, which may damage the lead. Use a new lead.

7. Verify lead placement with standard imaging techniques.
8. Coil excess lead in the subgaleal pocket (refer to Figure 11).

**Figure 11.** Coiling the lead in the subgaleal pocket.

9. If a second lead is implanted, repeat the target localization, lead implant, and test stimulation procedures. Tunnel the second lead based on the type of neurostimulator that will be implanted:
  - Tunnel the second lead on the opposite side of the skull from the first lead for bilateral stimulation with two single-program neurostimulators.
  - Tunnel the second lead next to the first lead for bilateral stimulation with a multi-program neurostimulator that will be implanted on one side of the body. (See Figure 12 for lead placement when a multi-program neurostimulator will be implanted.)

**Figure 12.** Lead tunneling when a multi-program neurostimulator is implanted.

**Note:** Place a ligature loosely on one of the leads to help identify the right or left brain lead position when two leads have been tunneled on the same side of the skull.

10. The remainder of the neurostimulation system may be implanted immediately after lead implant, or it may be implanted at a later date.
  - If implanted immediately, prepare the patient for implant of the extension and the neurostimulator.
  - If implanted at a later date, cap the lead(s) with the lead cap provided in the lead package.


**Note:** For instructions on implanting the extension and the neurostimulator, refer to the manuals for those components.

## Capping the lead

If the remainder of the neurostimulation system is not implanted immediately after lead implantation, perform the following steps:


1. Push the connector boot over the exposed end of the lead (Figure 13).


**Figure 13.** Push connector boot over lead.

2. Insert the exposed end of the lead into the lead cap contained in the lead kit (Figure 14).

**Figure 14.** Insert lead into lead cap.

3. Align the number 3 lead contact with the setscrew connector block (Figure 14).
4. Hold the setscrew connector block firmly between thumb and two fingers to prevent rotation (Figure 15).


**Figure 15.** Hold setscrew connector block between thumb and two fingers.

5. Tighten the single setscrew in the setscrew socket on the number 3 lead contact by turning it clockwise with the torque wrench provided (Figure 16). Tighten the setscrew until you hear multiple clicks.

**Cautions:**


- Use only the torque wrench provided in package or in compatible accessory kit. Use of an incorrect wrench may result in damage to the lead or an insufficient connection between the lead and extension, which would require surgical replacement or revision.
- Ensure the torque wrench is perpendicular to the setscrew surface when tightening the setscrew. Failure to hold the wrench perpendicular while tightening the setscrew may result in undertightening or false tightening of the setscrew.

**Figure 16.** Tighten setscrew.

6. Slide the connector boot over the lead cap (Figure 17).

**Note:** If it is difficult to position the boot, only sterile water may be used as a lubricant.


**Figure 17.** Slide connector boot over lead cap.

7. Place non-absorbable ligatures around both ends of the connector boot in the ligature channels (Figure 18).


**Cautions:**

- Do not overtighten ligatures because damage may occur to either the connector boot or the lead, which would require surgical replacement.
- Do not use polypropylene ligature material on silicone components. Polypropylene may damage the component, resulting in component failure and surgical replacement.

**Figure 18.** Ligature lead/lead cap.

8. Use blunt dissection to create a 50-mm diameter subgaleal pocket and a path to the desired lead-extension connection site (Figure 19).

**Figure 19.** Lead in subgaleal pocket.

9. Deliver the capped end of the lead to the desired temporary location (Figure 19).

**Cautions:**

- Be extremely careful when using sharp instruments around the lead body to avoid nicking or damaging the lead, which would require surgical replacement.
- Do not bend, kink, stretch, or twist the lead body when delivering it to the desired temporary location. Lead damage may occur, which would require surgical replacement.
- To avoid excessive traction when handling the lead with forceps, use only rubber-tipped bayonet forceps and grasp only the lead cap. Using the forceps on the lead body may cause lead damage and surgical replacement of the lead

10. In the subgaleal pocket, coil excess lead in a circle greater than 25 mm in diameter (Figure 19). Avoid bending or kinking the lead.
11. Close the incision and apply the appropriate dressing.

## Removing the lead cap

For each implanted lead:

1. Locate the lead cap at the proximal end of the lead.


**Note:** Use fluoroscopic observation or palpate the area with the fingers to locate the lead cap and confirm access without risk of cutting the lead.

2. Make an incision sufficient to expose the lead and lead cap. Allow room to hold the lead firmly when accessing the lead cap and when connecting the lead to the extension (Figure 20).

**Cautions:**


- Be extremely careful when using sharp instruments around the lead body to avoid nicking or damaging the lead, which would require surgical replacement.
- Do not pull the lead cap out of the tissue to prevent excess stress on the lead. If more slack is needed, gently grasp the lead body and lift the lead cap out of the tissue.
- Do not bend, kink, stretch, or twist the lead body when accessing the lead cap. Lead damage may occur, which would require surgical replacement.

**Figure 20.** Make incision to expose lead and lead cap.


3. Cut the distal ligature over the lead end of the connector boot (Figure 21).

△ **Caution:** Do not cut the lead when removing the ligature. Cutting the lead insulation can result in loss of stimulation and lead failure, which would require surgical replacement.


**Figure 21.** Cut distal ligature over lead end of connector boot.

4. Cut the connector boot over the lead cap to expose the setscrew (Figure 22).

**Figure 22.** Cut connector boot to expose setscrew.

5. Hold the setscrew connector block firmly between thumb and two fingers to prevent rotation (Figure 23).


**Figure 23.** Hold setscrew connector block between thumb and two fingers.

- Using the torque wrench, loosen the setscrew in the setscrew connector block by turning the wrench counterclockwise (approximately one turn) (Figure 24).

**Note:** Ensure the torque wrench is perpendicular to the setscrew surface while loosening the setscrew.


△ **Caution:** Hold the setscrew connector block firmly between thumb and two fingers to support the lead body and prevent rotation. Failure to properly prevent rotation may break or damage the lead, and may require lead replacement.

**Figure 24.** Loosen setscrew.


- Stabilize the lead and gently remove the lead from the setscrew connector block (Figure 25).

△ **Caution:** If resistance is felt while removing the lead from the lead cap, loosen the setscrew slightly to ensure that it clears the lead contact. Avoid disengaging the setscrew. Inspect the lead for damage (flattening or stretching of lead) if resistance was felt prior to removal, which would require surgical replacement.

**Figure 25.** Remove lead from setscrew connector block.

- Withdraw the lead cap and connector boot through the incision and discard (Figure 26).

**Figure 26.** Withdraw lead cap and connector boot.

**Tunneling to the neurostimulator site**

For instructions on tunneling to the neurostimulator site, refer to the appropriate extension manual.

**Making the lead-extension connection**

For instructions on making the lead-extension connection, refer to the appropriate extension manual.

**Making the extension-neurostimulator connection**

For instructions on making the extension-neurostimulator connection, refer to the appropriate neurostimulator manual.

**Performing intraoperative test stimulation with the lead-extension connection**

For instructions on making the lead-extension connection for intraoperative test stimulation, refer to the appropriate percutaneous extension manual.

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All Rights Reserved
Medtronic B.V.
Medtronic International Ltd.